



Anti-ZNF827 (aa 1027-1076) polyclonal antibody (CABT-BL3978)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 1027-1076 (VCKTKNMFER HLQIHLITRM FECDVCHKFM KTPEQLLEHK KCHTVPTGGL) of Human ZNF827.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [152485](#)

Protein Refseq [NP_849157](#)

UniProt ID [Q17R98](#)
