



Anti-ZNF835 (aa 108-157) polyclonal antibody (CABT-BL3980)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 108-157 (PGESPKERHP DSRQRE RGGG PKKPWKCGDC GKAFSYCSAF ILHQRIHTGE0 of Human ZNF835
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID [90485](#)

Protein Refseq [NP_001005850](#)

UniProt ID [Q9Y2P0](#)
