



Anti-ZNF860 (aa 608-632) polyclonal antibody (CABT-BL3983)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 608-632 of Human ZNF860.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	344787
Protein Refseq	NP_001131146

UniProt ID

[A6NHJ4](#)
