



Anti-ZNF879 polyclonal antibody (CABT-BL3985)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: MKSGGTNAGG SARETRRLSG AQAQLRPAAT GNVSELGVPL ASCAVWSCAP(Human)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	345462
Protein Refseq	NP_001129588

UniProt ID

[B4DU55](#)
