



Anti-ZNF880 (aa 51-81) polyclonal antibody (CABT-BL3986)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide between N terminal amino acids 51-81 of Human ZNF880, conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	400713
Protein Refseq	NP_001138906
UniProt ID	Q6PDB4