



Anti-ZNHIT3 (aa 25-75) polyclonal antibody (CABT-BL3990)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region between residues 25 and 75 of human ZNHIT3
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IP
Cellular Localization	Intracellular
Format	Liquid
Buffer	8mM PBS, 60mM Citrate, 150mM Tris, pH 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

GENE INFORMATION

Entrez Gene ID	9326
Protein Refseq	NP_004764

UniProt ID

[Q15649](#)
