



Anti-ZRANB3 (aa 1029-1079) polyclonal antibody (CABT-BL3996)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region between residue 1029 and 1079 of Human ZRANB3
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Buffer	Tris citrate/phosphate, pH 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

GENE INFORMATION

Entrez Gene ID	84083
Protein Refseq	NP_115519
UniProt ID	Q5FWF4