



## Anti-ZSCAN5D (aa 448-497) polyclonal antibody (CABT-BL3998)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: RKNLNEHKLI HSGEKPYKCP KCLRAFRPE TLKYHQKTHQ ETTAPRECEG, corresponding to amino acids 448-497 of Human LOC646698
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### GENE INFORMATION

Entrez Gene ID [646698](#)