



Anti-ZSWIM3 (aa 36-85) polyclonal antibody (CABT-BL3999)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 36-85 (VRFHNLNHGT SIREILYVQ VKFVCIRTQS NRKRTREADM CPAYLLLRN) of Human ZSWIM3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	140831
Protein Refseq	NP_542790

UniProt ID

[Q96MP5](#)
