



Anti-Postn polyclonal antibody (CABT-BL4000)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Mouse Periostin / POSTN protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	ELISA
Buffer	0.2 µm filtered solution in PBS with 5% trehalose
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -70°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

BACKGROUND

Introduction	Binds to heparin. Induces cell attachment and spreading and plays a role in cell adhesion. May play a role in extracellular matrix mineralization.
---------------------	--

GENE INFORMATION

Entrez Gene ID	50706
-----------------------	-----------------------

Protein Refseq [NP_001185694](#)

UniProt ID [Q96MP5](#)
