



Anti-C monoclonal antibody, clone 304633 (CABT-BL484)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	304633
Conjugate	Unconjugated
Applications	WB
Preservative	None
Storage	4°C or -20°C for longer storage

BACKGROUND

Introduction Beta-catenin is a transcriptional activator and oncoprotein involved in the development of several cancers. The protein encoded by this gene interacts directly with the C-terminal region of beta-catenin, inhibiting oncogenic beta-catenin-mediated transcriptional activation by competing with transcription factors for binding to beta-catenin. Two transcript variants encoding different isoforms have been found for this gene.

GENE INFORMATION

Entrez Gene ID [25776](#)

Protein Refseq [NP_001002880](#)

UniProt ID [A0A087WJW0](#)
