



Anti-RNF146 (full length) polyclonal antibody (CABT-BL5546)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse polyclonal antibody to Human RNF146.
Immunogen	Full length human RNF146 protein.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Sequence Similarities	Contains 1 RING-type zinc finger.Contains 1 WWE domain.
Format	Liquid
Size	50 µg
Buffer	PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.

GENE INFORMATION

Entrez Gene ID [81847](#)

Protein Refseq	NP_001229773
UniProt ID	Q9NTX7
Chromosome Location	6q22.1-q22.33
Function	ligase activity; metal ion binding; poly-ADP-D-ribose binding; protein binding; ubiquitin-protein ligase activity; zinc ion binding;