



Anti-ZC3H15 polyclonal antibody (CABT-BL5994)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-ZC3H15 polyclonal Antibody
Target	ZC3H15
Immunogen	Recombinant protein corresponding to amino acids of human ZC3H15.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Antigen affinity purification
Conjugate	Unconjugated
Sequence Similarities	LNELFKPVVAAQKISKGADPKSVVCAFFKQGQCTKGDCKKFSHDLTLERKCEKRSVYIDARDEELEKDTMDNWDEKKLEEVVNKKHGE
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	55854
Protein Refseq	NP_060941
UniProt ID	Q8WU90
Chromosome Location	2q32.1
Function	metal ion binding; nucleic acid binding; zinc ion binding;