



## Anti-Psph polyclonal antibody (CABT-BL6195)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Immunogen</b>	Antiserum is developed in rabbit using pure cultures of bacteria as immunogen.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Purification</b>	caprylic acid purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Size</b>	200 µl, 1 ml
<b>Buffer</b>	In 0.1 M phosphate buffered saline, pH 7.4.
<b>Preservative</b>	None
<b>Storage</b>	The antibody should be stored at 2-8°C. For storage longer than one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-life of 2 years after date of purchase.

### BACKGROUND

<b>Introduction</b>	<i>Pseudomonas syringae</i> pv. <i>phaseolicola</i> causes halo blight of the common bean, <i>Phaseolus vulgaris</i> , worldwide and remains difficult to control.
<b>Keywords</b>	<i>Pseudomonas syringae</i> pv <i>phaseolicola</i> ; Psph