



Anti-CERV polyclonal antibody [AP] (CABT-BL6227)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	AP
Applications	ELISA
Format	Liquid
Size	200 µl, 1 ml
Buffer	In 0.05 M TRIS buffer, pH 8.0.
Preservative	None
Storage	The antibody should be stored at 2-8°C. For storage longer than one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-life of 2 years after date of purchase.

BACKGROUND

Introduction	A DNA-containing virus with isometric particles about 45 nm in diameter occurring in cytoplasmic inclusion bodies. Transmissible by mechanical inoculation and by Myzus persicae. It is found naturally only in Dianthus caryophyllus. Distribution world-wide.
Keywords	Carnation etched ring virus; CERV