



# Anti-CVMV polyclonal antibody [AP] (CABT-BL6232)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	AP
Applications	ELISA
Format	Liquid
Size	200 µl, 1 ml
Buffer	In 0.05 M TRIS buffer, pH 8.0.
Preservative	None
Storage	The antibody should be stored at 2-8°C. For storage longer then one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-live of 2 years after date of purchase.

## BACKGROUND

Introduction	An RNA-containing virus with flexuous filamentous particles c. 790 nm long transmitted by the aphid Myzus persicae in the non-persistent manner, and by inoculation of sap. It causes disease in carnation and sweet william. Though uncommon in N.W. Europe, it is apparently more prevalent in USA and in some Mediterranean countries.
--------------	---

**Keywords**

Carnation vein mottle virus; CVMV

---