



Anti ABIN1 polyclonal antibody (CABT-BL6312)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	detects ABIN1 at ~73 kDa
Immunogen	mouse ABIN1 (residues 1-560)
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Mouse
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

BACKGROUND

Introduction This gene encodes an A20-binding protein which plays a role in autoimmunity and tissue homeostasis through the regulation of nuclear factor kappa-B activation. Mutations in this gene have been associated with psoriatic arthritis, rheumatoid arthritis, and systemic lupus erythematosus. Multiple transcript variants encoding different isoforms have been found for this

GENE INFORMATION

Entrez Gene ID [10318](#)

UniProt ID [B7Z8K2](#)
