



# Anti-UBXD7 polyclonal antibody (CABT-BL6316)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	detects UBXD7 at ~55 kDa
Immunogen	human UBXD7 (residues 1-489) [GST-tagged]
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

## GENE INFORMATION

Entrez Gene ID [26043](#)

UniProt ID

[O94888](#)

---