



# Anti-DCNL3 polyclonal antibody (CABT-BL6320)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                               |
|--------------------|-----------------------------------------------|
| Specificity        | : detects DCNL3 at<br>~34 kDa                 |
| Immunogen          | human DCNL3 (residues 1–105) [GST-tagged]     |
| Isotype            | IgG                                           |
| Source/Host        | Sheep                                         |
| Species Reactivity | Human                                         |
| Purification       | affinity-purified using immobilized immunogen |
| Conjugate          | Unconjugated                                  |
| Applications       | WB, IP                                        |
| Format             | Liquid                                        |
| Size               | 100 µg                                        |
| Buffer             | phosphate-buffered saline                     |
| Preservative       | None                                          |
| Storage            | 12 months at -20°C; aliquot as required       |

## GENE INFORMATION

Entrez Gene ID [123879](#)

UniProt ID

[Q8IWE4](#)

---