



Anti-DCNL5 polyclonal antibody (CABT-BL6323)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	detects DCNL5 at ~28 kDa
Immunogen	mouse DCNL5 (residue 1-237) [GST-tagged]
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Mouse
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

GENE INFORMATION

Entrez Gene ID [84259](#)

UniProt ID

[A0A024R3A2](#)
