



Anti-TRIAD1 polyclonal antibody (CABT-BL6330)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	detects TRIAD1 at ~58 kDa
Immunogen	human TRIAD1 (residues 1-493)
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

GENE INFORMATION

Entrez Gene ID [10425](#)

UniProt ID

[Q95376](#)
