



# Anti-TRIM65 polyclonal antibody (CABT-BL6331)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	detects TRIM65 at ~57 kDa
Immunogen	human TRIM65 (residues 1-485) [GST-tagged]
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

## GENE INFORMATION

Entrez Gene ID [201292](#)

UniProt ID

[Q6PJ69](#)

---