



# Anti-CYLD polyclonal antibody (CABT-BL6332)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	detects CYLD at ~110 kDa
Immunogen	human CYLD (residues 2-956) [His-tagged]
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

## BACKGROUND

Introduction	This gene encodes a cytoplasmic protein with three cytoskeletal-associated protein-glycine-conserved (CAP-GLY) domains that functions as a deubiquitinating enzyme. Mutations in this gene have been associated with cylindromatosis, multiple familial trichoepithelioma, and Brooke-Spiegler syndrome. Alternate transcriptional splice variants, encoding different
--------------	--

isoforms, have been characterized. [provided by RefSeq, Jul 2008]

---

## GENE INFORMATION

Entrez Gene ID [1540](#)

---

UniProt ID [A0A024R6T1](#)

---