



# Anti-OTULIN polyclonal antibody (CABT-BL6334)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	detects OTULIN at ~40 kDa
Immunogen	human OTULIN (residues 1-352)
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

## BACKGROUND

Introduction	0
--------------	---

# GENE INFORMATION

Entrez Gene ID [90268](#)

---

UniProt ID [Q96BN8](#)

---