



Anti-MAPK15 polyclonal antibody (CABT-BL6342)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------|---|
| Specificity | detects MAPK15 at ~63 kDa |
| Immunogen | human MAPK15 (residues 1-544) [6His-tagged] |
| Isotype | IgG |
| Source/Host | Sheep |
| Species Reactivity | Human |
| Purification | affinity-purified using immobilized immunogen |
| Conjugate | Unconjugated |
| Applications | WB |
| Format | Liquid |
| Size | 100 µg |
| Buffer | phosphate-buffered saline |
| Preservative | None |
| Storage | 12 months at -20°C; aliquot as required |

BACKGROUND

| | |
|--------------|-------------------------------------|
| Introduction | mitogen-activated protein kinase 15 |
|--------------|-------------------------------------|

GENE INFORMATION

Entrez Gene ID [225689](#)

UniProt ID [Q8TD08](#)
