



Anti-PINK1 (phospho T257 aa 250–262) polyclonal antibody (CABT-BL6344)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	detects PINK1 at ~63 kDa
Immunogen	human PINK1 (residues 250-262) [CAGEYGAV(pT)YRKSKR]
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

BACKGROUND

Introduction	This gene encodes a serine/threonine protein kinase that localizes to mitochondria. It is thought to protect cells from stress-induced mitochondrial dysfunction. Mutations in this gene cause one
--------------	--

GENE INFORMATION

Entrez Gene ID [65018](#)

UniProt ID [Q9BXM7](#)
