



Anti-C4ORF49 (N-terminal) polyclonal antibody (CABT-BL746)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from the N terminal domain of Human Ovary-specific acidic protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	IHC-P, WB
Cellular Localization	Membrane.
Format	Liquid
Size	200 µl
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [84709](#)

Protein Refseq [NP_116012](#)

UniProt ID [Q8TDB4](#)
