



Anti-C5ORF13 (C-terminal) polyclonal antibody (CABT-BL747)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from the C terminal residues of Human P311.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Cellular Localization	Cytoplasmic
Format	Liquid
Buffer	1% BSA, 0.1M Tris glycine, pH 7.2
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Neuronal protein 3.1 is a protein that in humans is encoded by the C5orf13 gene.

GENE INFORMATION

Entrez Gene ID [9315](#)

Protein Refseq [NP_001135953](#)

UniProt ID [Q16612](#)
