



## Anti-C5ORF55 (aa 2-51) polyclonal antibody (CABT-BL754)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 2-51 (PAVFMLASSS ALQCGRGVPR FPRTEVGAGH SVNEETKAEK VGNQTSVIPA) of human C5orf55.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Secreted
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

### GENE INFORMATION

Entrez Gene ID	<a href="#">116349</a>
Protein Refseq	<a href="#">NP_612473</a>