



# Anti-C6ORF10 (aa 513-562) polyclonal antibody (CABT-BL756)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide corresponding to a region within C terminal amino acids 513-562 (KKEKGDKN TK GDKGKDKVKG KRESEINGEK SKGSKRAKAN TGRKYNKKVE) of Human C6ORF10.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Membrane.
<b>Format</b>	Liquid
<b>Buffer</b>	2% Sucrose, PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction** Uncharacterized protein C6orf10 is a protein that in humans is encoded by the C6orf10 gene.

# GENE INFORMATION

Entrez Gene ID [10665](#)

---

Protein Refseq [NP\\_006772](#)

---

UniProt ID [Q13647](#)

---