



# Anti-IL17A Monoclonal antibody, clone 75EFD28 [FITC] (CABT-BL7718)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	75EFD28
Purification	Affinity chromatography
Conjugate	FITC
Applications	FC
Format	Liquid
Concentration	5ul/Test
Size	25 tests
Buffer	PBS, pH 7.2, with 0.1% gelatin, 0.2% BSA
Preservative	0.09% Sodium Azide
Storage	4°C, store in dark, DO NOT FREEZE!

## BACKGROUND

Introduction	Interleukin-17A (IL-17A) is a CD4+ T cell-derived cytokine that promotes inflammatory responses in cell lines and is elevated in rheumatoid arthritis, asthma, multiple sclerosis,
--------------	--

psoriasis, and transplant rejection. The cDNA encoding human IL-17A was isolated from a library of CD4+ T cells; the encoded protein exhibits 72 percent amino acid identity with HVS13 , an open reading frame from a T lymphotropic Herpesvirus saimiri, and 63 percent with mouse CTLA-8 (cytotoxic T-lymphocyte associated antigen-8). Human IL-17A exists as glycosylated 20-30 kD homodimers. High levels of IL-17A homodimer are produced by activated peripheral blood CD4+ T-cells. IL-17A enhances expression of the intracellular adhesion molecule-1 (ICAM-1) in human fibroblasts. Human IL-17A also stimulates epithelial, endothelial, or fibroblastic cells to secrete IL-6, IL-8, G-CSF, and PGE2. In the presence of human IL-17A, fibroblasts can sustain the proliferation of CD34+ hematopoietic progenitors and induce maturation into neutrophils. Mouse, rat, and human IL-17A can induce IL-6 secretion in mouse stromal cells, indicating that all homologs can recognize the mouse receptor. IL-17A binds to a receptor that binds also to HVS13 (viral IL-17A) and to CTLA-8.

---

<b>Keywords</b>	IL17A; interleukin 17A; interleukin-17A; cytotoxic T-lymphocyte-associated antigen 8; IL17; Ctl8; IL-17; Ctl8; IL-17A;
-----------------	--

---

## GENE INFORMATION

---

<b>Entrez Gene ID</b>	<a href="#">3605</a>
-----------------------	----------------------

---

<b>UniProt ID</b>	<a href="#">Q16552</a>
-------------------	------------------------

---