



# Anti-C9ORF150 (aa 11-40) polyclonal antibody (CABT-BL807)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, between amino acids 11-40 from the N-terminal region of Human C9orf150.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

## GENE INFORMATION

Entrez Gene ID	<a href="#">286343</a>
Protein Refseq	<a href="#">NP_981948</a>
UniProt ID	<a href="#">Q8IV03</a>