



Anti-C9ORF30 (full length) polyclonal antibody (CABT-BL810)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length Human C9orf30 protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Rat, Human
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The function of this protein remains unknown.
--------------	---

GENE INFORMATION

Entrez Gene ID	91283
----------------	-----------------------

Protein Refseq [NP_001185734](#)

UniProt ID [A0A024R171](#)
