



Anti-C9ORF4 (internal region) polyclonal antibody (CABT-BL811)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal amino acids 215-264, (HDDNGRVRIQHFYNVGQWAKEIQRNPARDEEGVFENNRV TCRFKRPVNVP) of Human C9ORF4
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [23732](#)

Protein Refseq [NP_055149](#)

UniProt ID [Q9P0K9](#)
