



Anti-C9ORF64 (aa 266-315) polyclonal antibody (CABT-BL813)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to C-terminal amino acids 266-315 (EMLSYGDRQE VEIRGCSLWC VELIRDCLLE LIEQKGEKPN GEINSILLDY) of Human C9orf64.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	98% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	84267
Protein Refseq	NP_115683
UniProt ID	Q5T6V5