



# Anti-C9ORF68 (full length) polyclonal antibody (CABT-BL814)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Full length Human C9orf68 protein (AAH34293).
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	1X PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

Entrez Gene ID	<a href="#">55064</a>
Protein Refseq	<a href="#">NP_001034484</a>
UniProt ID	<a href="#">Q8N4H0</a>