



Anti-C9ORF68 (aa 190-382) polyclonal antibody (CABT-BL815)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment, corresponding to a region within the internal sequence amino acids 190-382 of Human C9orf68 (AAH34293)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	20% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	55064
Protein Refseq	NP_001034484
UniProt ID	Q8N4H0