



Anti-C9ORF78 (aa 1-209) polyclonal antibody (CABT-BL816)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment corresponding to a region within amino acids 1 - 209 of Human C9orf78.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC/IF
Format	Liquid
Size	50 µl
Buffer	10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction Uncharacterized protein C9orf78 is a protein that in humans is encoded by the C9orf78 gene.

GENE INFORMATION

Entrez Gene ID [51759](#)

Protein Refseq [NP_057604](#)

UniProt ID [Q9NZ63](#)
