



Anti-CACNG8 (aa 180-229) polyclonal antibody (CABT-BL822)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 180-229 (ISANAGEPGP KRDEEKKNHY SYGWSFYFGG LSFILAEVIG VLAVNIYIER) of Mouse CACNG8.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane. Cell junction; synapse; postsynaptic cell membrane; postsynaptic density.
Format	Liquid
Buffer	97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [81905](#)

Protein Refseq [NP_573453](#)

UniProt ID [Q8VHW2](#)
