



Anti-CAF1-105 (aa 610-709) polyclonal antibody (CABT-BL830)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment, corresponding to amino acids 610-709 of Fruit fly (<i>Drosophila melanogaster</i>) p105 CAF1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, ICC/IF, IP
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	36107
Protein Refseq	NP_610589
UniProt ID	A1Z898