



Anti-CAF1-180 polyclonal antibody (CABT-BL831)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment (Fruit fly (Drosophila melanogaster))
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	31801
Protein Refseq	NP_572495
UniProt ID	Q9W3D1