



# Anti-CALY polyclonal antibody (CABT-BL833)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                                                                                     |
|--------------------|-----------------------------------------------------------------------------------------------------|
| Immunogen          | Recombinant protein corresponding to amino acids of human CALY.                                     |
| Isotype            | IgG                                                                                                 |
| Source/Host        | Rabbit                                                                                              |
| Species Reactivity | Human                                                                                               |
| Purification       | Antigen affinity purification                                                                       |
| Conjugate          | Unconjugated                                                                                        |
| Applications       | IHC                                                                                                 |
| Format             | Liquid                                                                                              |
| Buffer             | In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)                                                   |
| Preservative       | 0.02% Sodium Azide                                                                                  |
| Storage            | Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing. |

## BACKGROUND

|              |                                                                                                                                     |
|--------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Introduction | Neuron-specific vesicular protein calcyon is a protein that in humans is encoded by the CALY gene. Its alternative name is Calcyon. |
|--------------|-------------------------------------------------------------------------------------------------------------------------------------|

## GENE INFORMATION

Entrez Gene ID [50632](#)

---

Protein Refseq [NP\\_056537](#)

---

UniProt ID [Q9NYX4](#)

---