



# Anti-CAMTA1 polyclonal antibody (CABT-BL835)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human CAMTA1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC, IF
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	CAMTA1 is a transcription activator which is a putative tumor suppressor. It is thought to interact with calmodulin. Defects in CAMTA1 are detected in patients with oligodendroglioma and astrocytoma.
--------------	---

# GENE INFORMATION

Entrez Gene ID [23261](#)

---

Protein Refseq [NP\\_001182492](#)

---

UniProt ID [Q9Y6Y1](#)

---