



# Anti-CAPS polyclonal antibody (CABT-BL843)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Full length calcyphosine protein (Human)
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasm.
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

Introduction	This gene encodes a calcium-binding protein, which may play a role in the regulation of ion transport. A similar protein was first described as a potentially important regulatory protein in the dog thyroid and was termed as R2D5 antigen in rabbit. Alternative splicing of this gene generates two transcript variants.
--------------	--

# GENE INFORMATION

Entrez Gene ID	<a href="#">828</a>
Protein Refseq	<a href="#">NP_004049</a>
UniProt ID	<a href="#">Q13938</a>