



Anti-CAR1 (aa 2-51) polyclonal antibody (CABT-BL847)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the N terminal amino acids 2-51 (ASPDWGYDDKNGPEQWSKLYPIANGNNQSPVDIKTSETK HDTSLKPISVS) of Human Carbonic Anhydrase I.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasm.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [12346](#)

Protein Refseq [NP_001077426](#)

UniProt ID [P13634](#)
