



Anti-CASKIN1 (aa 1416-1430) polyclonal antibody (CABT-BL849)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide (aa 1416 - 1430 of CASKIN 1) coupled to key-hole limpet hemocyanin via an added N-terminal cysteine residue.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, ICC
Preservative	None
Storage	-20°C, Avoid freeze / thaw cycles

GENE INFORMATION

Entrez Gene ID	57524
Protein Refseq	NP_065815.1
UniProt ID	Q8WXD9