



Anti-CCDC101 (aa 188-218) polyclonal antibody (CABT-BL854)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 188-218 of Human SGF29, conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Cellular Localization	Nuclear
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	CCDC101 is a subunit of 2 histone acetyltransferase complexes: the ADA2A (TADA2A; MIM 602276)-containing (ATAC) complex; and the SPT3 (SUPT3H; MIM 602947)-TAF9 (MIM 600822)-GCN5 (KAT2A; MIM 602301)/PCAF (KAT2B; MIM 602303) acetylase; (STAGA)
---------------------	---

complex. Both of these complexes contain either GCN5 or PCAF, which are paralogous acetyltransferases (Wang et;al., 2008 (PubMed 18838386)).

GENE INFORMATION

Entrez Gene ID	112869
----------------	------------------------

Protein Refseq	NP_612423
----------------	---------------------------

UniProt ID	Q96ES7
------------	------------------------
