



Anti-CCDC165 (aa 225-275) polyclonal antibody (CABT-BL859)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region between amino acids 225-275 of Human KIAA0802.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IP
Format	Liquid
Buffer	Tris citrate/phosphate, pH 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	KIAA0802 also known as hypothetical protein LOC23255 is a protein that is encoded by the KIAA0802 gene.
--------------	---

GENE INFORMATION

Entrez Gene ID [23255](#)

Protein Refseq [NP_056025](#)

UniProt ID [Q9Y4B5](#)
