



Anti-CCDC62 (aa 654-684) polyclonal antibody (CABT-BL865)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 654-684 of Human CCDC62, conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasm. Nucleus. Note: Mainly nuclear
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	84660
Protein Refseq	NP_958843

UniProt ID

[Q6P9F0](#)
