



Anti-CCDC67 (aa 36-85) polyclonal antibody (CABT-BL867)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within N terminal amino acids 36-85 (WERKMRALETRLDLRDQELANAQTCLDQKGQEVGLLRQK LDSLEKCNLAM) of Human CCDC67, NP_857596
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	159989
Protein Refseq	NP_857596

UniProt ID

[Q05D60](#)
